

## LAMTOR5 抗原（重组蛋白）

中文名称：LAMTOR5 抗原（重组蛋白）

英文名称：LAMTOR5 Antigen (Recombinant Protein)

别名：XIP; HBXIP

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	late endosomal/lysosomal adaptor, MAPK and MTOR activator 5
<b>Synonyms:</b>	XIP; HBXIP
<b>Swissprot:</b>	O43504
<b>Gene Accession:</b>	BC062619
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that specifically complexes with the C-terminus of hepatitis B virus X protein (HBx). The function of this protein is to negatively regulate HBx activity and thus to alter the replication life cycle of the virus.