

兔抗 BCL2(Phospho-Ser70) 多克隆抗体

中文名称：兔抗 BCL2(Phospho-Ser70) 多克隆抗体

英文名称： Anti-BCL2(Phospho-Ser70) rabbit polyclonal antibody

别名： Bcl-2; PPP1R50

储存： 冷冻（-20℃） 避光

抗原： BCL2(Phospho-Ser70)

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Two transcript variants, produced by alternate splicing, differ in their C-terminal ends.
Applications:	IHC
Name of antibody:	BCL2(Phospho-Ser70)
Immunogen:	Peptide sequence around phosphorylation site of serine 70(R-T-S(p)-P-L) derived from Human BCL-2.

Full name:	B-cell CLL/lymphoma 2
Synonyms :	Bcl-2; PPP1R50
SwissProt:	P10415
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

