

兔抗 BCL2(Phospho-Ser87) 多克隆抗体

中文名称：兔抗 BCL2(Phospho-Ser87) 多克隆抗体

英文名称： Anti-BCL2(Phospho-Ser87) rabbit polyclonal antibody

别名： Bcl-2; PPP1R50

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： BCL2(Phospho-Ser87)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Two transcript variants, produced by alternate splicing, differ in their C-terminal ends.
Applications:	WB
Name of antibody:	BCL2(Phospho-Ser87)
Immunogen:	Peptide sequence around phosphorylation site of serine 87 (A-L-S(p)-P-V) derived from Human BCL-2.

Full name:	B-cell CLL/lymphoma 2
Synonyms :	Bcl-2; PPP1R50
SwissProt:	P10415
WB Predicted band size:	26 kDa
WB Positive control:	Hela cells lysates treated with Nocodazole
WB Recommended dilution:	500-1000

