

LASP1 抗原（重组蛋白）

中文名称：LASP1 抗原（重组蛋白）

英文名称：LASP1 Antigen (Recombinant Protein)

别名：MLN50; Lasp-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	LIM and SH3 protein 1
Synonyms:	MLN50; Lasp-1
Swissprot:	Q14847
Gene Accession:	BC012460
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a subfamily of LIM proteins, characterized by a LIM motif and a domain of Src homology region 3, and also a member of the nebulin family of actin-binding proteins. The encoded protein is a cAMP and cGMP dependent signaling protein and binds to the actin cytoskeleton at extensions of the cell membrane. The encoded protein has been linked to metastatic breast cancer, hematopoietic tumors such as B-cell lymphomas, and colorectal cancer.