

LGALS13 抗原（重组蛋白）

中文名称：LGALS13 抗原（重组蛋白）

英文名称： LGALS13 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	lectin, galactoside-binding, soluble, 13
Synonyms:	PP13, GAL13, PLAC8
Swissprot:	Q9UHV8
Gene Accession:	BC069312
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. The protein encoded by this gene has lysophospholipase activity. It is composed of two identical subunits which are held together by disulfide bonds. This protein has structural similarity to several members of the beta-galactoside-binding S-type lectin family.