

## 兔抗 BCL2L14 多克隆抗体

中文名称：兔抗 BCL2L14 多克隆抗体

英文名称： Anti-BCL2L14 rabbit polyclonal antibody

别名： BCLG

相关类别： 一抗

抗原： BCL2L14

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BCL2L14
<b>Immunogen:</b>	Fusion protein of human BCL2L14

<b>Full name:</b>	BCL2-like 14 (apoptosis facilitator)
<b>Synonyms :</b>	BCLG
<b>SwissProt:</b>	Q9BZR8
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and human tonsil
<b>IHC Recommend dilution:</b>	50-200

