

LYPD4 抗原（重组蛋白）

中文名称：LYPD4 抗原（重组蛋白）

英文名称：LYPD4 Antigen (Recombinant Protein)

别名：SMR

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	LY6/PLAUR domain containing 4
Synonyms:	SMR
Swissprot:	Q6UWN0
Gene Accession:	BC034629
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The specific function of LYPD4 is not yet known.