

兔抗 CuNVgp5 多克隆抗体

- 中文名称：兔抗 CuNVgp5 多克隆抗体
- 英文名称：Anti-CuNVgp5 rabbit polyclonal antibody
- 别名：p20 protein
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：CuNVgp5
- 反应种属：Cucumber necrosis virus
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	plays a role in initiation of gene silencing at the early stage of viral infection.
Applications:	WB
Name of antibody:	CuNVgp5
Immunogen:	Synthetic peptide of Cucumber necrosis virus CuNVgp5
Full name:	CuNVgp5
Synonyms :	p20 protein
SwissProt:	P15184