

兔抗 BCL2 多克隆抗体

中文名称：兔抗 BCL2 多克隆抗体

英文名称：Anti-BCL2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：BCL2

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Two transcript variants, produced by alternate splicing, differ in their C-terminal ends.
Applications:	ELISA, WB, IHC
Name of antibody:	BCL2
Immunogen:	Synthetic peptide of human BCL2
Full name:	B-cell CLL/lymphoma 2
SwissProt:	P10415
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human colorectal cancer and Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	22 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-2000



