

## 兔抗 CXADR 多克隆抗体

中文名称：兔抗 CXADR 多克隆抗体

英文名称： Anti-CXADR rabbit polyclonal antibody

别名： CAR; HCAR; CAR4/6

抗原： CXADR

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CXADR
<b>Immunogen:</b>	Synthesized peptide derived from internal of human CXADR.
<b>Full name:</b>	coxsackie virus and adenovirus receptor
<b>Synonyms :</b>	CAR; HCAR; CAR4/6

<b>SwissProt:</b>	P78310
<b>WB Predicted band size:</b>	40 kDa
<b>WB Positive control:</b>	COLO cells and HepG2 cells lysates
<b>WB Recommended dilution:</b>	500-3000

