

## MAPK12 抗原（重组蛋白）

中文名称：MAPK12 抗原（重组蛋白）

英文名称：MAPK12 Antigen (Recombinant Protein)

别名：ERK3; ERK6; SAPK3; PRKM12; SAPK-3; P38GAMMA

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	mitogen-activated protein kinase 12
<b>Synonyms:</b>	ERK3; ERK6; SAPK3; PRKM12; SAPK-3; P38GAMMA
<b>Swissprot:</b>	P53778
<b>Gene Accession:</b>	BC015741
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Activation of members of the mitogen-activated protein kinase family is a major mechanism for transduction of extracellular signals. Stress-activated protein kinases are one subclass of MAP kinases. The protein encoded by this gene functions as a signal transducer during differentiation of myoblasts to myotubes.