

## 兔抗 CXCL3 多克隆抗体

中文名称：兔抗 CXCL3 多克隆抗体

英文名称： Anti-CXCL3 rabbit polyclonal antibody

别名： GRO3; GROg; MIP2B; SCYB3; MIP-2b; CINC-2b

抗原： CXCL3

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Chemokine (C-X-C motif) ligand 3 (CXCL3) is a small cytokine belonging to the CXC chemokine family that is also known as GRO3 oncogene (GRO3), GRO protein gamma a (GROg) and macrophage inflammatory protein-2-beta (MIP2b). CXCL3 controls migration and adhesion of monocytes and mediates its effects on its target cell by interacting with a cell surface chemokine receptor called CXCR 2. More recently, it has been shown that Cxcl3 regulates cell autonomously the migration of the precursors of cerebellar granule neurons toward the internal layers of cerebellum, during the morphogenesis of cerebellum.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CXCL3
<b>Immunogen:</b>	Fusion protein of human CXCL3
<b>Full name:</b>	chemokine (C-X-C motif) ligand 3
<b>Synonyms :</b>	GRO3; GROg; MIP2B; SCYB3; MIP-2b; CINC-2b
<b>SwissProt:</b>	P19876
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

