

兔抗 BCL9 多克隆抗体

中文名称: 兔抗 BCL9 多克隆抗体

英文名称: Anti-BCL9 rabbit polyclonal antibody

别名: LGS

储存: 冷冻 (-20℃)

抗原: BCL9

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	BCL9 is associated with B-cell acute lymphoblastic leukemia. It may be a target of translocation in B-cell malignancies with abnormalities of 1q21. Its function is unknown. The overexpression of BCL9 may be of pathogenic significance in B-cell malignancies.
Applications:	ELISA, IHC
Name of antibody:	BCL9
Immunogen:	Fusion protein of human BCL9
Full name:	B-cell CLL/lymphoma 9
Synonyms:	LGS
SwissProt:	O00512

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	25-100



