

## MED18 抗原（重组蛋白）

中文名称：MED18 抗原（重组蛋白）

英文名称：MED18 Antigen (Recombinant Protein)

别名：p28b

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	mediator complex subunit 18
<b>Synonyms:</b>	p28b
<b>Swissprot:</b>	Q9BUE0
<b>Gene Accession:</b>	BC002694
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).