

MIF4GD 抗原（重组蛋白）

中文名称：MIF4GD 抗原（重组蛋白）

英文名称：MIF4GD Antigen (Recombinant Protein)

别名：MIFD; AD023; SLIP1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	MIF4G domain containing
Synonyms:	MIFD; AD023; SLIP1
Swissprot:	A9UHW6
Gene Accession:	BC033759
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein which interacts with the N-terminus of the stem-loop binding protein (SLBP) and the 3' end of histone mRNA. This interaction facilitates the activation of histone mRNA translation. Alternative splicing results in multiple transcript variants encoding different isoforms.