

兔抗 CXCR5 多克隆抗体

中文名称：兔抗 CXCR5 多克隆抗体

英文名称： Anti-CXCR5 rabbit polyclonal antibody

别名： C-X-C motif chemokine receptor 5; BLR1; CD185; MDR15

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CXCR5

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a multi-pass membrane protein that belongs to the CXC chemokine receptor family. It is expressed in mature B-cells and Burkitt's lymphoma. This cytokine receptor binds to B-lymphocyte chemoattractant (BLC), and is involved in B-cell migration into B-cell follicles of spleen and Peyer patches. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CXCR5
Immunogen:	Synthetic peptide of human CXCR5

Full name:	C-X-C motif chemokine receptor 5
Synonyms:	BLR1; CD185; MDR15
SwissProt:	P32302
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human colorectal cancer
IHC Recommend dilution:	50-200



