

MOB1B 抗原（重组蛋白）

中文名称：MOB1B 抗原（重组蛋白）

英文名称：MOB1B Antigen (Recombinant Protein)

别名：MATS2; MOB4A; MOBKL1A

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	MOB kinase activator 1B
Synonyms:	MATS2; MOB4A; MOBKL1A
Swissprot:	Q7L9L4
Gene Accession:	BC038112
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is similar to the yeast Mob1 protein . Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Three transcript variants encoding different isoforms have been found for this gene.