

## MOB4 抗原（重组蛋白）

中文名称： MOB4 抗原（重组蛋白）

英文名称： MOB4 Antigen (Recombinant Protein)

别名： 2C4D; MOB1; MOB3; PHOCN; PREI3; CGI-95; MOBKL3

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	MOB family member 4, phocein
<b>Synonyms:</b>	2C4D; MOB1; MOB3; PHOCN; PREI3; CGI-95; MOBKL3
<b>Swissprot:</b>	Q9Y3A3
<b>Gene Accession:</b>	BC005237
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene was identified based on its similarity with the mouse counterpart. Studies of the mouse counterpart suggest that the expression of this gene may be regulated during oocyte maturation and preimplantation following zygotic gene activation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HSPE1.