

MRPL1 抗原(重组蛋白)

中文名称: MRPL1 抗原(重组蛋白)

英文名称: MRPL1 Antigen (Recombinant Protein)

别 名: mitochondrial ribosomal protein L1; L1MT; BM022; MRP-L1

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 51-325 amino acids of human MRPL1

技术规格

Full name:	mitochondrial ribosomal protein L1
Synonyms:	L1MT; BM022; MRP-L1
Swissprot:	Q9BYD6
Gene Accession:	BC014356
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitocho ndrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA comp osition compared to prokaryotic ribosomes, where this ratio is reverse d. Another difference between mammalian mitoribosomes and prokary otic ribosomes is that the latter contain a 5S rRNA. Among different s



pecies, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L1 ribosomal protein family.