

兔抗 CYC1 多克隆抗体

中文名称：兔抗 CYC1 多克隆抗体

英文名称：Anti-CYC1 rabbit polyclonal antibody

别名：UQCR4; MC3DN6

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CYC1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a subunit of the cytochrome bc1 complex, which plays an important role in the mitochondrial respiratory chain by transferring electrons from the Rieske iron-sulfur protein to cytochrome c. Mutations in this gene may cause mitochondrial complex III deficiency, nuclear type 6.
Applications:	WB
Name of antibody:	CYC1
Immunogen:	Synthesized peptide derived from internal of human CYC1.
Full name:	cytochrome c-1
Synonyms :	UQCR4; MC3DN6

SwissProt:	P08574
WB Predicted band size:	35 kDa
WB Positive control:	HepG2 cells lysate
WB Recommended dilution:	500-3000

