

## 兔抗 BHLHA9 多克隆抗体

中文名称：兔抗 BHLHA9 多克隆抗体

英文名称： Anti-BHLHA9 rabbit polyclonal antibody

别名： basic helix-loop-helix family member a9; BHLHF42

相关类别： 一抗

抗原： BHLHA9

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the basic helix-loop-helix family . The encoded protein is a transcription factor involved in limb development. Mutations in this gene have been associated with mesoaxial synostotic syndactyly Malik-Percin type (MSSD). Copy number variation of a locus containing this gene has been linked to a form of split-hand/foot malformation with long bone deficiency (SHFLD3).
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	BHLHA9
<b>Immunogen:</b>	Synthetic peptide of human BHLHA9
<b>Full name:</b>	basic helix-loop-helix family member a9

<b>Synonyms:</b>	BHLHF42
<b>SwissProt:</b>	Q7RTU4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	Human cerebrum tissue lysate
<b>WB Recommended dilution:</b>	200-1000

