

## MYL6 抗原（重组蛋白）

中文名称：MYL6 抗原（重组蛋白）

英文名称：MYL6 Antigen (Recombinant Protein)

别名：LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	myosin, light chain 6, alkali, smooth muscle and non-muscle
<b>Synonyms:</b>	LC17; ESMLC; LC17A; LC17B; MLC-3; MLC1SM; MLC3NM; MLC3SM; LC17-GI; LC17-NM
<b>Swissprot:</b>	P60660
<b>Gene Accession:</b>	BC017455
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain that is expressed in smooth muscle and non-muscle tissues. Genomic sequences representing several pse

udogenes have been described and two transcript variants encoding different isoforms have been identified for this gene.