

## MYOZ1 抗原（重组蛋白）

中文名称：MYOZ1 抗原（重组蛋白）

英文名称：MYOZ1 Antigen (Recombinant Protein)

别名：CS-2, FATZ, MYOZ

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	myozenin 1
<b>Synonyms:</b>	CS-2, FATZ, MYOZ
<b>Swissprot:</b>	Q9NP98
<b>Gene Accession:</b>	BC025753
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is primarily expressed in the skeletal muscle, and belongs to the myozenin family. Members of this family function as calcineurin-interacting proteins that help tether calcineurin to the sarcomere of cardiac and skeletal muscle. They play an important role in modulation of calcineurin signaling.