

兔抗 BICD1 多克隆抗体

中文名称：兔抗 BICD1 多克隆抗体

英文名称：Anti-BICD1 rabbit polyclonal antibody

别名：BICD

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BICD1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an adaptor protein that belongs to the bicaudal D family of dynein cargo adaptors. The encoded protein acts as an intracellular cargo transport cofactor that regulates the microtubule-based loading of cargo onto the dynein motor complex. It also controls dynein motor activity and coordination. It has a domain architecture consisting of coiled-coil domains at the N- and C-termini that are highly conserved in other family members. Naturally occurring mutations in this gene are associated with short telomere length and emphysema.
Applications:	ELISA, IHC

Name of antibody:	BICD1
Immunogen:	Synthetic peptide of human BICD1
Full name:	BICD cargo adaptor 1
Synonyms:	BICD; bic-D 1
SwissProt:	Q96G01
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

