

## 兔抗 CYLD 多克隆抗体

中文名称：兔抗 CYLD 多克隆抗体

英文名称：Anti-CYLD rabbit polyclonal antibody

别名：EAC; MFT; SBS; TEM; BRSS; CDMT; MFT1; CYLD1; CYLDI; USPL2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CYLD

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cytoplasmic protein with three cytoskeletal-associated protein-glycine-conserved (CAP-GLY) domains that functions as a deubiquitinating enzyme. Mutations in this gene have been associated with cylindromatosis, multiple familial trichoepithelioma, and Brooke-Spiegler syndrome. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CYLD
<b>Immunogen:</b>	Synthetic peptide of human CYLD
<b>Full name:</b>	cylindromatosis (turban tumor syndrome)

<b>Synonyms:</b>	EAC; MFT; SBS; TEM; BRSS; CDMT; MFT1; CYLD1; CYLDI; USPL2
<b>SwissProt:</b>	Q9NQC7
<b>IHC positive control:</b>	Human gastric cancer and human brain
<b>IHC Recommend dilution:</b>	50-200

