

兔抗 CYP1A1/2 多克隆抗体

中文名称：兔抗 CYP1A1/2 多克隆抗体

英文名称： Anti-CYP1A1/2 rabbit polyclonal antibody

别名： AHH; AHRR; CP11; CYP1; P1-450; P450-C; P450DX/CP12; P3-450; P450(PA)

抗原： CYP1A1/2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acids, and xenobiotics.
Applications:	WB
Name of antibody:	CYP1A1/2
Immunogen:	Synthesized peptide derived from internal of human Cytochrome P450 1A1/2.
Full name:	cytochrome P450, family 1, subfamily A, polypeptide 1/2
Synonyms :	AHH; AHRR; CP11; CYP1; P1-450; P450-C; P450DX/CP12; P3-

	450; P450(PA)
SwissProt:	P04798/P05177
WB Predicted band size:	58 kDa
WB Positive control:	COS-7 cells, RAW264.7 cells and HT-29 cells lysates
WB Recommended dilution:	500-3000

