

兔抗 BICRA 多克隆抗体

中文名称：兔抗 BICRA 多克隆抗体

英文名称：Anti-BICRA rabbit polyclonal antibody

别名：BRD4 interacting chromatin remodeling complex associated protein; GLTSCR1

储存：冷冻 (-20℃)

抗原：BICRA

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in BRD4-mediated gene transcription.
Applications:	ELISA, IHC
Name of antibody:	BICRA
Immunogen:	Synthetic peptide of human BICRA
Full name:	BRD4 interacting chromatin remodeling complex associated protein
Synonyms:	GLTSCR1
SwissProt:	Q9NZM4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	40-200

