

## 兔抗 BIN3 多克隆抗体

中文名称：兔抗 BIN3 多克隆抗体

英文名称： Anti-BIN3 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： BIN3

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene is a member of the BAR domain protein family. The encoded protein is comprised solely of a BAR domain which is predicted to form coiled-coil structures and proposed to mediate dimerization, sense and induce membrane curvature, and bind small GTPases. BAR domain proteins have been implicated in endocytosis, intracellular transport, and a diverse set of other processes
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	BIN3
<b>Immunogen:</b>	Fusion protein of human BIN3
<b>Full name:</b>	bridging integrator 3
<b>SwissProt:</b>	Q9NQY0
<b>ELISA Recommended dilution</b>	1000-2000

:	
<b>WB Predicted band size:</b>	30 kDa
<b>WB Positive control:</b>	231, A549, LoVo and HepG2 cells
<b>WB Recommended dilution:</b>	200-1000

