

NDUFB10 抗原（重组蛋白）

中文名称：NDUFB10 抗原（重组蛋白）

英文名称：NDUFB10 Antigen (Recombinant Protein)

别名：PDSW

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 10, 22kDa
Synonyms:	PDSW
Swissprot:	O96000
Gene Accession:	BC007509
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 10 is an enzyme that in humans is encoded by theNDUFB10 gene. NADH dehydrogenase (ubiquinone) 1 beta subcomplex subunit 10 is an accessory subunit of the NADH dehydrogenase (ubiquinone) complex, located in the mitochondrial inner membrane. It is also known as Complex I and is the largest of the five complexes of the electron transport chain.