

## 兔抗 CYP2B6 多克隆抗体

- 中文名称：兔抗 CYP2B6 多克隆抗体
- 英文名称：Anti-CYP2B6 rabbit polyclonal antibody
- 别名：CPB6; EFVM; IIB1; P450; CYP2B; CYP2B7; CYP2B7P; CYP11B6
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：CYP2B6
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene, CYP2B6, encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and its expression is induced by phenobarbital. The enzyme is known to metabolize some xenobiotics, such as the anti-cancer drugs cyclophosphamide and ifosfamide.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CYP2B6

<b>Immunogen:</b>	Synthetic peptide of human CYP2B6
<b>Full name:</b>	cytochrome P450, family 2, subfamily B, polypeptide 6
<b>Synonyms:</b>	CPB6; EFVM; IIB1; P450; CYP2B; CYP2B7; CYP2B7P; CYP2B6
<b>SwissProt:</b>	P20813
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	56 kDa
<b>WB Positive control:</b>	Human colon cancer tissue
<b>WB Recommended dilution:</b>	200-1000



