

兔抗 BLK 多克隆抗体

中文名称：兔抗 BLK 多克隆抗体

英文名称： Anti-BLK rabbit polyclonal antibody

别名： MODY11

相关类别： 一抗

抗原： BLK

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a nonreceptor tyrosine-kinase of the src family of proto-oncogenes that are typically involved in cell proliferation and differentiation. The protein has a role in B-cell receptor signaling and B-cell development. The protein also stimulates insulin synthesis and secretion in response to glucose and enhances the expression of several pancreatic beta-cell transcription factors.
Applications:	ELISA, WB, IHC
Name of antibody:	BLK
Immunogen:	Fusion protein of human BLK
Full name:	B lymphoid tyrosine kinase

Synonyms :	MODY11
SwissProt:	P51451
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	58 kDa
WB Positive control:	A172 cells
WB Recommended dilution:	200-1000



