

兔抗 BLNK(Phospho-Tyr96) 多克隆抗体

中文名称：兔抗 BLNK(Phospho-Tyr96) 多克隆抗体

英文名称： Anti-BLNK(Phospho-Tyr96) rabbit polyclonal antibody

别名： bca; AGM4; BASH; LY57; SLP65; BLNK-S; SLP-65

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： BLNK(Phospho-Tyr96)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a cytoplasmic linker or adaptor protein that plays a critical role in B cell development. This protein bridges B cell receptor-associated kinase activation with downstream signaling pathways, thereby affecting various biological functions. The phosphorylation of five tyrosine residues is necessary for this protein to nucleate distinct signaling effectors following B cell receptor activation.
Applications:	WB
Name of antibody:	BLNK(Phospho-Tyr96)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 96(D-S-Y(p)-E-P) derived from Human BLNK .

Full name:	B-cell linker
Synonyms :	bca; AGM4; BASH; LY57; SLP65; BLNK-S; SLP-65
SwissProt:	Q8WV28
WB Predicted band size:	50 kDa
WB Positive control:	COLO205 lysates
WB Recommended dilution:	500-1000

