

兔抗 CYP2U1 多克隆抗体

中文名称：兔抗 CYP2U1 多克隆抗体

英文名称： Anti-CYP2U1 rabbit polyclonal antibody

别名： SPG49; SPG56; P450TEC

储存： 冷冻（-20℃） 避光

抗原： CYP2U1

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This enzyme is a hydroxylase that metabolizes arachidonic acid, docosahexaenoic acid, and other long chain fatty acids.
Applications:	WB
Name of antibody:	CYP2U1
Immunogen:	Synthesized peptide derived from internal of human Cytochrome P450 2U1.

Full name:	cytochrome P450, family 2, subfamily U, polypeptide 1
Synonyms :	SPG49; SPG56; P450TEC
SwissProt:	Q7Z449
WB Predicted band size:	62 kDa
WB Positive control:	Hela cells and Lovo cells lysates
WB Recommended dilution:	500-3000

