

## 兔抗 CYP2W1 多克隆抗体

中文名称：兔抗 CYP2W1 多克隆抗体

英文名称： Anti-CYP2W1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： CYP2W1

宿 主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CYP2W1
<b>Immunogen:</b>	Fusion protein of human CYP2W1
<b>Full name:</b>	cytochrome P450, family 2, subfamily W, polypeptide 1
<b>SwissProt:</b>	Q8TAV3
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human cervical cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	100-300

