

兔抗 BLOC1S6 多克隆抗体

中文名称：兔抗 BLOC1S6 多克隆抗体

英文名称： Anti-BLOC1S6 rabbit polyclonal antibody

别名： biogenesis of lysosomal organelles complex 1 subunit 6; PA; HPS9; PLDN; BLOS6; PALLID

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： BLOC1S6

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene may play a role in intracellular vesicle trafficking. It interacts with Syntaxin 13 which mediates intracellular membrane fusion. Mutations in this gene cause symptoms associated with Hermansky-Pudlak syndrome-9. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.
Applications:	ELISA, WB, IHC
Name of antibody:	BLOC1S6
Immunogen:	Fusion protein of human BLOC1S6

Full name:	biogenesis of lysosomal organelles complex 1 subunit 6
Synonyms:	PA; HPS9; PLDN; BLOS6; PALLID
SwissProt:	Q9UL45
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	20 kDa
WB Positive control:	231, HeLa, Jurkat, 293T, K562 and HepG2 cell lysates
WB Recommended dilution:	1000-5000



