

NUDT21 抗原（重组蛋白）

中文名称：NUDT21 抗原（重组蛋白）

英文名称：NUDT21 Antigen (Recombinant Protein)

别名：CPSF5; CFIM25

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	nudix (nucleoside diphosphate linked moiety X)-type motif 21
Synonyms:	CPSF5; CFIM25
Swissprot:	O43809
Gene Accession:	BC001403
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. This gene encodes the 25kD subunit of the protein complex, which is composed of four polypeptides.