

ORMDL3 抗原（重组蛋白）

中文名称：ORMDL3 抗原（重组蛋白）

英文名称：ORMDL3 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

| | |
|---------------------------|--|
| Full name: | ORMDL sphingolipid biosynthesis regulator 3 |
| Swissprot: | Q8N138 |
| Gene Accession: | BC000638 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | Negative regulator of sphingolipid synthesis. May indirectly regulate endoplasmic reticulum-mediated Ca(+2) signaling. SNPs on 17q21 locus that are associated with childhood asthma also show a consistent and strong association with transcript levels of ORMDL3, indicating that genetic variants regulating ORMDL3 expression are determinants of susceptibility to childhood asthma. |