

兔抗 CYP4A11/22 多克隆抗体

中文名称：兔抗 CYP4A11/22 多克隆抗体

英文名称：Anti-CYP4A11/22 rabbit polyclonal antibody

别名：CP4Y; CYP4A2; CYP4A11

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CYP4A11/22

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and hydroxylates medium-chain fatty acids such as laurate and myristate.
Applications:	WB
Name of antibody:	CYP4A11/22
Immunogen:	Synthesized peptide derived from internal of human CYTOCHROME P450 4A11/22.

Full name:	cytochrome P450, family 4, subfamily A, polypeptide 11/22
Synonyms :	CP4Y; CYP4A2; CYP4AII
SwissProt:	Q02928/Q5TCH4
WB Predicted band size:	59 kDa
WB Positive control:	COLO205, Huvec and Hela cells lysates
WB Recommended dilution:	500-3000

