

## OTULINL 抗原（重组蛋白）

中文名称：OTULINL 抗原（重组蛋白）

英文名称：OTULINL Antigen (Recombinant Protein)

别名：OTU deubiquitinase with linear linkage specificity like; NET20; FAM105A

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 157-356 amino acids of human OTULINL

技术规格

<b>Full name:</b>	OTU deubiquitinase with linear linkage specificity like
<b>Synonyms:</b>	NET20; FAM105A
<b>Swissprot:</b>	Q9NUU6
<b>Gene Accession:</b>	BC011524
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Lacks deubiquitinase activity.