

兔抗 CYP4F2 多克隆抗体

中文名称：兔抗 CYP4F2 多克隆抗体

英文名称： Anti-CYP4F2 rabbit polyclonal antibody

别名： CPF2

抗原： CYP4F2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	is gene encodes a member of the cytochrome P450 super family of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum. The enzyme starts the process of inactivating and degrading leukotriene B4, a potent mediator of inflammation. This gene is part of a cluster of cytochrome P450 genes on chromosome 19. Another member of this family, CYP4F11, is approximately 16 kb away.
Applications:	IHC

Name of antibody:	CYP4F2
Immunogen:	Synthesized peptide derived from N-terminal of human Cytochrome P450 4F2.
Full name:	cytochrome P450, family 4, subfamily F, polypeptide 2
Synonyms :	CPF2
SwissProt:	P78329
IHC positive control:	Human liver carcinoma tissue
IHC Recommend dilution:	50-100

