

兔抗 BMP8A 多克隆抗体

中文名称：兔抗 BMP8A 多克隆抗体

英文名称： Anti-BMP8A rabbit polyclonal antibody

抗 原： BMP8A

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Induces cartilage and bone formation. May be the osteoinductive factor responsible for the phenomenon of epithelial osteogenesis. Plays a role in calcium regulation and bone homeostasis. By similarity. Signaling protein involved in regulation of thermogenesis and energy balance. Proposed to increase the peripheral response of brown adipose tissue (BAT) to adrenergic stimulation while acting centrally in the hypothalamus to increase sympathetic output to BAT.
Applications:	WB, IF
Name of antibody:	BMP8A
Immunogen:	Synthesized peptide derived from C-terminal of human BMP8A.

Full name:	bone morphogenetic protein 8a
Synonyms :	0
SwissProt:	Q7Z5Y6
WB Predicted band size:	45 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000
IF positive control:	A549 cells
IF Recommend dilution:	100-500



