

## 兔抗 CYP8B1 多克隆抗体

中文名称：兔抗 CYP8B1 多克隆抗体

英文名称： Anti-CYP8B1 rabbit polyclonal antibody

别名： CP8B; CYP12

抗原： CYP8B1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	Involved in bile acid synthesis and is responsible for the conversion of 7 alpha-hydroxy-4-cholesten-3-one into 7 alpha, 12 alpha-dihydroxy-4-cholesten-3-one. Responsible for the balance between formation of cholic acid and chenodeoxycholic acid. Has a rather broad substrate specificity including a number of 7-alpha-hydroxylated C27 steroids.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CYP8B1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human Cytochrome P450 8B1.
<b>Full name:</b>	cytochrome P450, family 8, subfamily B, polypeptide 1

<b>Synonyms :</b>	CP8B; CYP12
<b>SwissProt:</b>	Q9UNU6
<b>IHC positive control:</b>	Human liver carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000



