

兔抗 BMP10 多克隆抗体

中文名称：兔抗 BMP10 多克隆抗体

英文名称： Anti-BMP10 rabbit polyclonal antibody

抗 原： BMP10

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse, Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the TG F-beta family of growth factors. Data suggest that the similar protein in mouse plays an important role in trabeculation of the embryonic heart. In human, this protein may signal through receptor serine/threonine kinases.
Applications:	ELISA, WB, IHC
Name of antibody:	BMP10
Immunogen:	Synthetic peptide of human BMP10
Full name:	bone morphogenetic protein 10
SwissProt:	O95393
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

WB Predicted band size:	48 kDa
WB Positive control:	Human breast cancer tissue lysate
WB Recommended dilution:	500-2000



