

## PDLIM4 抗原（重组蛋白）

中文名称：PDLIM4 抗原（重组蛋白）

英文名称：PDLIM4 Antigen (Recombinant Protein)

别名：RIL

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	PDZ and LIM domain 4
<b>Synonyms:</b>	RIL
<b>Swissprot:</b>	P50479
<b>Gene Accession:</b>	BC016765
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis.