

PHB 抗原（重组蛋白）

中文名称：PHB 抗原（重组蛋白）

英文名称：PHB Antigen (Recombinant Protein)

别名：PHB1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	prohibitin
Synonyms:	PHB1
Swissprot:	P35232
Gene Accession:	BC013401
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is evolutionarily conserved, and its product is proposed to play a role in human cellular senescence and tumor suppression. Antiproliferative activity is reported to be localized to the 3' UTR, which is proposed to function as a trans-acting regulatory RNA. Several pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants