

兔抗 BNIP2 多克隆抗体

中文名称：兔抗 BNIP2 多克隆抗体

英文名称：Anti-BNIP2 rabbit polyclonal antibody

别名：NIP2; BNIP-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BNIP2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector.
Applications:	ELISA, IHC
Name of antibody:	BNIP2
Immunogen:	Synthetic peptide of human BNIP2
Full name:	BCL2/adenovirus E1B 19kDa interacting protein 2
Synonyms :	NIP2; BNIP-2

SwissProt:	Q12982
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and Human gastric cancer
IHC Recommend dilution:	50-200

