

兔抗 BNIP3L 多克隆抗体

中文名称：兔抗 BNIP3L 多克隆抗体

英文名称：Anti-BNIP3L rabbit polyclonal antibody

别名：NIX; BNIP3a

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BNIP3L

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. This protein counteracts the apoptotic inducer BNIP3 and may play a role in tumor suppression.
Applications:	ELISA, IHC
Name of antibody:	BNIP3L
Immunogen:	Synthetic peptide of human BNIP3L

Full name:	BCL2/adenovirus E1B 19kDa interacting protein 3-like
Synonyms:	NIX; BNIP3a
SwissProt:	O60238
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	15-50

