

兔抗 CYP26A1 多克隆抗体

中文名称：兔抗 CYP26A1 多克隆抗体

英文名称： Anti-CYP26A1 rabbit polyclonal antibody

别名： CP26; CYP26; P450RAI; P450RA1

抗原： CYP26A1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Plays a key role in retinoic acid metabolism. Acts on retinoids, including all-trans-retinoic acid (RA) and its stereo isomer 9-cis-RA. Capable of both 4-hydroxylation and 18-hydroxylation. Responsible for generation of several hydroxylated forms of RA, including 4-OH-RA, 4-oxo-RA and 18-OH-RA.
Applications:	WB
Name of antibody:	CYP26A1
Immunogen:	Synthesized peptide derived from internal of human Cytochrome P450 26A1.
Full name:	cytochrome P450, family 26, subfamily A, polypeptide 1

Synonyms :	CP26; CYP26; P450RAI; P450RAI1
SwissProt:	O43174
WB Predicted band size:	56 kDa
WB Positive control:	HepG2 cells, Hela cells and 293 cells lysates
WB Recommended dilution:	500-3000

