

PRSS38 抗原（重组蛋白）

中文名称：PRSS38 抗原（重组蛋白）

英文名称：PRSS38 Antigen (Recombinant Protein)

别名：MPN2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	protease, serine, 38
Synonyms:	MPN2
Swissprot:	A1L453
Gene Accession:	BC130400
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein remains unknown.