

## PRTFDC1 抗原（重组蛋白）

中文名称：PRTFDC1 抗原（重组蛋白）

英文名称： PRTFDC1 Antigen (Recombinant Protein)

别 名： HHGP

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	phosphoribosyl transferase domain containing 1
<b>Synonyms:</b>	HHGP
<b>Swissprot:</b>	Q9NRG1
<b>Gene Accession:</b>	BC008662
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	PRTFDC1 (phosphoribosyltransferase domain-containing protein 1), also known as HHGP, is a 225 amino acid protein that belongs to the purine/pyrimidine phosphoribosyltransferase family. The functions of PRTFDC1 are not well understood; however it has been suggested that it is frequently silenced in oral squamous-cell carcinomas by aberrant promoter hypermethylation. Conserved in vertebrates, PRTFDC1 is suggested to be a possible tumor-suppressor. PRTFDC

1 is encoded by a gene which maps to chromosome 10.